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The following is an addendum to the **THREE DIMENSIONAL HOFFMAN BRAIN PHANTOM User Manual**. This information supersedes descriptions of “cover plate”, “cover plate removal screws”, **Parts List**, and **assembly and disassembly procedures**.

## HOFFMAN

### **Assembly Procedure:** Hoffman 3D Brain, new style

General suggestions: Always start all screws in a hole pattern by hand before tightening any with a screwdriver. Thumb screws can be tightened firmly by hand, but never with tools such as pliers. Be sure there is a light film of silicone lubricant around each O-ring. Tighten screws having O-rings gently and slowly to achieve a complete seal.

#### **I. 3D Brain Insert**

Parts required:

- nineteen Brain Slices (individual discs)
- acrylic Top Spacer Plate (nine holes)
- acrylic Bottom Spacer Plate (fifteen holes)
- four 6” Delrin through-bolts

Procedure

1. Stack the slices so the scribed lines on their edges make an uninterrupted “V” pattern. This is a guide for correct sequencing.
2. Locate the hole in the Top Spacer Plate that is nearest the scribed lines, and insert a through-bolt with its head in the counterbored hole. Rest the plate on a clean table, bolt pointed upward.
3. One by one, gently slide the slices onto the bolt, maintaining continuity of the scribed lines.
4. Slide the assembly carefully out over the edge of the table, and thread the bolt loosely into the acrylic bottom plate, finishing the “V” pattern.
5. Rotate the assembly 180 degrees, insert a bolt in the opposite hole, and wiggle it through all of the slices. Thread it into the Bottom plate.
6. Carefully flip the assembly and insert the remaining two bolts. Slowly and evenly tighten the bolts with a large screwdriver, working around the pattern several times.

The 3D Brain Insert is now ready to be placed in the Cylinder.

To disassemble, remove all four through-bolts.

Important: Observe the following sequence of main assembly. These steps are all detailed in Section II below.

- A. Install bottom Center Plug in Cylinder
- B. Place Insert in Cylinder and tighten Bottom Screws
- C. Place Lid on Cylinder and install Top Screws

## **II. Main Assembly**

Parts Required:

- 3D Brain Insert
- Cylinder
- Center Bottom Plug with O-ring (large flanged head)
- four Slotted Bottom Screws with O-rings (blue)
- Lid with O-ring
- six Top Thumb Screws with O-rings
- five Filler Plugs with O-rings

Procedure

1. Thread the Center Bottom Plug into the cylinder and tighten gently with a large screwdriver or coin.
2. Place the cylinder on its side and carefully slide the Insert into it, Top Plate first. Note that top of head corresponds to *bottom* of Cylinder.
3. Rotate the Insert about the center plug to align the bolt patterns.
4. Thread-in the four Slotted Bottom Screws (blue) and tighten them gently with a large screwdriver.
5. Set the cylinder on its bottom. Position the Lid so that the bolt pattern corresponds to the pattern on Insert, then press it firmly onto the Cylinder with both hands. You may need to rotate the Lid so that the bolt patterns line up precisely.
6. Insert three of the large Filler Plugs in holes not needed for filling or venting, and screw each down until its O-ring contacts the Lid face.
7. Thread-in the six Top Thumb Screws. Tighten them slowly and firmly by hand, working around the pattern several times.

You are now ready to fill the Hoffman 3D Brain Phantom

Disassembly after draining is the reverse of this procedure; Remove Filler Plugs and Top Screws and take off Lid. Place Cylinder on side and remove Bottom Screws; then slide out Insert.